

## Trapeze Bar

The trapeze bar allows a patient more independence by providing the patient with a means of self help to change position in bed, move onto and off of a bedpan, or transfer into or out of a bed. A representative from our company will install the trapeze bar for you, insuring for tightness of clamps and stability.



### Warnings

- If any part of the trapeze bar becomes loose or unstable discontinue use and call our office. Continued use of the trapeze bar may result in serious injury.
- Check all clamps daily to ensure that they are securely tight.
- Adjust the trapeze bar height by hooking the chain on the "S" hook at the correct height.
- To get the bar out of the way simply loop it over the horizontal bar of the trapeze frame.

### Warranty

If you purchased your trapeze bar the warranty is for 1 year from initial set up.

If you are renting your trapeze bar the warranty is for the length of the rental period.